goldtrace

Goldtrace Fetal Spiral Electrode (FSE)

REF CNS000004

EN	Instructions for Use	2
FR	Notice d'utilisation	4
DE	Gebrauchsanweisung	6
NL	Gebruiksaanwijzing	8
ES	Instrucciones de uso	10
PT	Instruções de utilização	12
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	16
SV	Bruksanvisning	18
FI	Käyttöohjeet	20
DA	Brugsanvisning	22
HU	Használati utasítás	24
CS	Návod k použití	26
SK	Návod na použitie	28
RO	Instrucțiunile de utilizare	30
ET	Kasutusjuhend	32
LT	Naudojimo instrukcija	34
LV	Lietošanas pamācība	36
EL	Οδηγίες χρήσης	38
PL	Instrukcja użytkowania	40
HR	Upute za uporabu	42
BG	Инструкции за употреба	44
SL	Navodila za uporabo	46
IS	Notkunarleiðbeiningar	48

EN Instructions for Use

Goldtrace Fetal Spiral Electrode (FSE)

Device description

Goldtrace FSE is a sterile, single use device transmitting the fetal ECG signal for continuous electronic intrapartum monitoring.

Intended purpose

The use of Goldtrace FSE is indicated for any labor and delivery requiring continuous electronic intrapartum monitoring of the fetal heart rate or analysis of the fetal ECG waveform. The Goldtrace FSE is intended to be used with a compatible adapter cable (Legplate) connector specified by Neoventa Medical AB.

Indications for use

Goldtrace FSE is indicated for patients requiring continuous electronic intrapartum monitoring of the fetal heart rate during labor and delivery.

Contraindications

Do not apply FSE to fetal face, fontanels, or genitalia, or when placenta previa is present or suspected, or when it is not possible to identify fetal presenting part where application is being considered. Do not apply FSE in the presence of active herpes lesions, hepatitis C or HIV infection or when woman is a confirmed carrier of hemophilia and fetus is affected or status is unknown. Use caution and refer to your own institution's policies and procedures when fetus is premature, or woman is positive for group B streptococcus, hepatitis A, hepatitis B, syphilis, or gonorrhea. FSE shall not be used during or in combination with operative delivery. In case operative delivery is needed, remove FSE.

Warnings

This FSE is not defibrillation-proof and should not be used with defibrillation-proof ECG connections.

Due to design of FSE, penetration of fetal epidermis may possibly cause trauma, hemorrhage and/or infection. FSE must be used under conditions of aseptic technique. Amniotic membranes must be ruptured prior to attachment of FSE. Do not use FSE as tool to rupture amniotic membranes. Remove FSE from patient before performing any electro-surgical procedures. Single use only. Reuse can cause contamination and infection. Do not use if packaging or sterile barrier is damaged before use. Do not pull the Spiral Needle from the fetal skin. Do not pull the Wires apart. Do not over-rotate the Spiral Needle during attachment.


Caution

Federal (US) law restricts this device to sale by or on the order of a physician.

Electrode-tissue impedance can potentially lead to degraded ECG signal quality. Users are advised to monitor signal quality and adjust usage conditions, as necessary.

Intended users

FSE is intended to be used by medically trained personnel at delivery wards.

Compatible devices

FECG Legplate Adapter for Neoventa

FECG Legplate Adapter for Philips

Philips DECG Reusable Legplate Adapter cable

Please refer to the Instructions for Use for the Legplate adapter cable and fetal monitoring system for information on reprocessing instructions. Please refer to the Instructions for Use of the fetal monitoring system for EMC related information.


For explanation of symbols used see *Symbols glossary*: <https://www.neoventa.com/symbols-glossary/>

Intended patient population


The intended patient population is the fetus during delivery and the mother.

Applying the FSE


1. Using aseptic technique, open the package and remove FSE. Leave the Wires locked in the Retention Notch and check that the Protection Tab is in locked position. Straighten the Wires if necessary.
2. Form the Guide Tube gently to the desired anatomical curvature.
3. Perform a vaginal examination, to clearly identify the fetal presenting part.
4. Insert FSE at correct angle to the fetal presenting part until the Guide Tube reaches the presenting part. Place the Guide Tube firmly against fetal presenting part.


5. To release the Protection Tab, retract the Drive Tube Grip away from the Guide Tube until the Protection Tab clears the Guide Tube (approximately 1 inch or 2.5 cm).


6. Maintain pressure against fetal presenting part with the Guide Tube. Advance the Drive Tube until the Spiral Needle reaches the presenting part.


7. With pressure against fetal presenting part, rotate the Drive Tube Grip clockwise until you feel full attachment. Note: "full attachment" usually occurs when you feel mild resistance to further rotation and recoil from the Drive Tube Grip. This typically happens after one full rotation. **DO NOT OVER-ROTATE.**

8. Remove the Wires from the Retention Notch. Slide the Drive Tube off the Wires. Then slide the Guide Tube off the Wires and wipe Connector clean.


9. Insert the Connector end securely into compatible adapter cable (Legplate) connector.


10. Once all preparations are completed, ensure signal quality.

Removing the FSE

1. Pull the Connector end of FSE out of the compatible adapter cable (Legplate) connector.
2. Remove the Spiral Needle by grasping the Hub and the Wires as close to fetal presenting part as possible and twist counter-clockwise until free from presenting part.

DO NOT PULL THE SPIRAL NEEDLE FROM THE FETAL SKIN. DO NOT PULL THE WIRES APART.

3. Inspect the Spiral Needle to ensure FSE is intact. If the Spiral Needle has separated from the Hub, remove it by using aseptic technique.

Disposal

The Spiral Needle shall be handled as sharp object. Follow institution's policies and procedures for safe disposal.

Reporting responsibilities

Serious incidents, related to FSE, shall be reported to manufacturer and competent authorities without undue delay.

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