

# goldtrace

## The Gold-plated Fetal Scalp Electrode

### A message from the heart

Goldtrace® Fetal Spiral Electrode (FSE) is used when fetal heart rate monitoring during labor requires continuous and reliable internal registration. Goldtrace FSE is made of materials whose properties ensure the best-possible-quality signal for CTG.

Goldtrace FSE comes in a flexible protective guide tube that keeps the spiral tip inside the tube during application, protecting the user, mother and baby. The flexibility of the guide tube makes it possible to form according to desired anatomical curvature.

### Specifications

- Gold-plated spiral electrode
- Single use only
- Includes disposable applicator
- Not made with natural rubber latex
- Sterile product



### Article number

CNS000004/75      Goldtrace Fetal Spiral Electrode  
in boxes of 75 single packed pcs

### Materials

|                 |                                       |
|-----------------|---------------------------------------|
| Spiral needle   | 316 Stainless Steel with 10-25µm Gold |
| Reference plate | 316 Stainless Steel with 10-25µm Gold |
| Hub             | PE                                    |
| Connector       | PE and brass, gold-plated             |
| Wires           | Copper with PVC jacket                |
| Drive tube      | PE                                    |
| Guide tube      | PE                                    |
| Drive tube grip | PE                                    |

Follow the link for more  
information: [https:// www.  
neovolta.com/products/goldtrace/](https://www.neovolta.com/products/goldtrace/)



PRM 400008/44 R5B © 2023 Neoventa Medical AB. All rights reserved.

#### Neoventa Medical AB

Norra Ågatan 32, 431 35 Mölndal, Sweden  
+46 31 758 32 00, [info@neoventa.com](mailto:info@neoventa.com)

[www.neoventa.com](http://www.neoventa.com)

**neoventa**   
A better start in life