

# Mounting Solutions

## Compatible with Stan S41 Maternal and Fetal Monitor

There are two mounting options for the S41 monitor: on a trolley or on a wall arm. Both options offer a flexible, practical, and ergonomic mounting solution. All wall arms and trolleys include an adaptor plate compatible with the S41 monitor and enables the quick-release function.

Pick and mix among trolleys, wall arms, and accessories to build a mounting system tailored entirely after your needs.

### Trolley

The trolley enables easy access to necessary equipment in different locations when needed, and meets various ergonomic needs with adjustable height and position options.

- Available in fixed or variable height
- Quick release function
- Enables free positioning of the monitor
- Stable with quiet-running, easy-rolling lockable wheels and anti-static wire
- Friction-resistant tilt functionality
- A variety of accessories available



#### Trolley Fixed Height

Trolley with fixed height and quick-release function for the S41 monitor.

- Height: 128 cm
- Maximum storage weight: 4,5 kg
- Color: Cream white

#### NVT-0002-60

- 1 Trolley fixed height
- 1 Bin
- 1 Handle
- 2 Cable management clips
- 1 Anti-static wire



#### Trolley Variable Height

Trolley with variable height (gas-spring assisted) and quick-release function for the S41 monitor.

- Height range: 120-156 cm
- Maximum storage weight: 6,8 kg
- Color: Cream white

#### SBS120430

- 1 Trolley variable height
- 1 Bin
- 1 Handle
- 2 Cable management clips
- 1 Anti-static wire



# Wall Arm

The wall arm offers increased free floor space, and stores accessories in a space-preserving and practical way for different ergonomic needs and room layouts.

- Available in Standard and Premium
- Quick release function
- Lateral and vertical movement functionality and tilt function
- Available in long and short wall channel lengths
- A variety of accessories available

Standard



Premium

---

## Standard Wall Arm – Short Channel

Fixed height wall arm with lateral movement and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 39 cm
- Wall channel length: 48,3 cm

## SBS120415/1

- 1 Wall arm lateral movement
- 1 Wall channel 48,3 cm



---

## Standard Wall Arm – Long Channel

Fixed height wall arm with lateral movement and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 39 cm
- Wall channel length: 63,5 cm

## SBS120415/2

- 1 Wall arm lateral movement
- 1 Wall channel 63,5 cm



---

## Premium Wall Arm – Short Channel

Wall arm with lateral and vertical movement (gas-spring assisted) and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 66 cm
- Wall channel length: 48,3 cm

## SBS120420/1

- 1 Wall arm lateral and vertical movement
- 1 Wall channel 48,3 cm



---

### Premium Wall Arm – Long Channel

Wall arm with lateral and vertical movement (gas-spring assisted) and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 66 cm
- Wall channel length: 63,5 cm

### SBS120420/2

- 1 Wall arm lateral and vertical movement
- 1 Wall channel 63,5 cm



---

## Accessories

Trolley: ● Wall arm: ●

---

### Basket

Storage for monitoring accessories such as cables, blood pressure cuffs.

Dimensions: W:32 cm D:19 cm H:15 cm

- RS-0001-24
- WM-0001-18B



---

### Cord wrap

Cord management to be mounted horizontally or vertically.

- RS-0013-03E



---

### Chart Bin

Chart bin with quick release adapter.

Dimension: W:36,2 cm. D:8,3 cm, H:22,9 cm

- WM-0001-23
- WM-0001-21



---

## Bin

● WM-0001-22

Bin with removable dividers. Quick release adapter for easy removal. Included in the trolley standard package.

● WM-0001-27

Dimensions: H:7,6-12,7 cm D:30,5 cm W:30 cm



---

## Quad 4-hook

● UT-0001-11

Hook for cable management.

Maximum load per hook: 1,27 kg each.



Follow link to find out more.

[www.neoventa.com](http://www.neoventa.com)



PRM 800401/EN R1C © 2024 Neoventa Medical AB. All rights reserved.

### Neoventa Medical AB

Norra Ågatan 32, 431 35 Mölndal, Sweden

+46 31 758 32 00, [info@neoventa.com](mailto:info@neoventa.com)

[www.neoventa.com](http://www.neoventa.com)

**neoventa**   
A better start in life