

# Accessories and Consumables

## Stan S41

---

### Accessories for Internal Fetal Monitoring

#### **FECG Legplate Adapter for Neoventa**

**ACC300110**

Adapter cable connecting a Fetal Spiral Electrode (FSE) and Skin electrode to a Fetal monitor for continuous electronic intrapartum monitoring.

Validated with Goldtrace Fetal Spiral Electrode (FSE).



#### **Goldtrace Fetal Spiral Electrode**

**CNS000004/75**

Fetal Spiral Electrode for internal fetal heart rate and FECG registration.

Single use, 75 pcs/unit



#### **Skin Electrode for FECG**

**CNS000003/30**  
**CNS000003/600**

Wet gel skin electrode from ConMed with characteristics suitable for use with fetal ST Analysis.

Single use.

Available in boxes of 30 and 600.



#### **IUP Adapter Cable for Utah Medical/Intran**

**P1263-03027**

For internal monitoring of uterine activity. Compatible with IntranPlus IUP catheters from Utah Medical.

#### **Abrasive Paper**

**CNS000401/1500**

Abrasive paper for skin preparation prior to application of skin electrode.

Single use, 1500 pcs/unit

# Accessories for External Fetal Monitoring

---

## Wired Fetal Movement Marker

P1221-12035

Remote with button for manual registration of fetal movements.



## Wireless Fetal Movement Marker

P1271-12006

Wireless remote with button for manual registration of fetal movements.



# Accessories for Maternal Monitoring

---

## 3-Leadwire ECG Cable

P9001-00478

Leadset for monitoring of ECG, heart rate, and respiratory rate.



## Skin Electrode for MECG Recording

P7001-00296

Skin electrodes for monitoring of MECG.

25 pcs/package.



## MSpO2 Sensor

P7002-00008

Sensor for registration of pulse and oxygen saturation.

Requires P9001-00501 "MSpO2 Extension Cable" for use.



## MSpO2 Extension Cable

P9001-00501

Extension cable for SpO2 with connector.

Requires P7002-00008 "MSpO2 fingerclip" for use.



## NIBP Cuff Extension Hose

P9001-00485

NIBP extension cable with connector.



---

**NIBP Cuff Adult****P9001-00473**

Cuff for measurement of blood pressure (perimeter 25-35 cm).



---

**NIBP Cuff Large Adult****P9001-00474**

Cuff for measurement of blood pressure (perimeter 33-47 cm).



---

**NIBP Cuff Adult Thigh****P9001-00482**

Cuff for measurement of blood pressure (perimeter 46-66 cm).



---

## Consumables

---

**Legplate Belt****CNS000106**

Stretch belt for the legplate.

- 3 cm wide
- Roll of 15 m
- Not made with natural rubber latex



---

**Transducer Belt – Regular****CNS000107**

Stretch belt for the TOCO and FHR transducers.

- 5 cm wide
- Roll of 15 m
- Not made with natural rubber latex



---

**Transducer Belt – Wide****CNS000108**

Stretch belt for the TOCO and FHR transducers.

- 10 cm wide
- Roll of 15 m
- Not made with natural rubber latex



---

**Transducer Belt – Precut****P2224-08001**

Reusable stretch belt for the TOCO and FHR transducers.

- 6 cm wide
- 150 cm long
- Pack of 2 pcs
- Not made with natural rubber latex




---

### Printer Paper (CTG and ST Analysis)

**P8105-00063**

Printer paper with CTG+ST analysis grid, with 50-210 bpm FHR scale and 20 bpm per cm scaling.

3 pcs/unit.



---

### Printer Paper (CTG Only)

**P8105-00003**

Printer paper with CTG-only grid, with 50-210 bpm FHR scale and 20 bpm per cm scaling.

3 pcs/unit.

