

# Accessories and Consumables

## Stan S41

---

### Accessories for Internal Fetal Monitoring

#### FECG Legplate Adapter for Neoventa

ACC300110

Adapter cable connecting a Fetal Spiral Electrode (FSE) and Skin electrode to a Fetal monitor for continuous electronic intrapartum monitoring.

Validated with Goldtrace Fetal Spiral Electrode (FSE).



#### Goldtrace Fetal Spiral Electrode

CNS000004/75

Fetal Spiral Electrode for internal fetal heart rate and FECG registration.

Single use, 75 pcs/unit



#### Skin Electrode for FECG

CNS000003/30  
CNS000003/600

Wet gel skin electrode from ConMed with characteristics suitable for use with fetal ST Analysis.

Single use.

Available in boxes of 30 and 600.



#### IUP Adapter Cable for Utah Medical/Intran

P1263-03027

For internal monitoring of uterine activity. Compatible with IntranPlus IUP catheters from Utah Medical.

#### Abrasive Paper

CNS000401/1500

Abrasive paper for skin preparation prior to application of skin electrode.

Single use, 1500 pcs/unit

# Accessories for External Fetal Monitoring

---

## Wired Fetal Movement Marker

P1221-12035

Remote with button for manual registration of fetal movements.



# Accessories for Maternal Monitoring

---

## 3-Leadwire ECG Cable

P9001-00478

Leadset for monitoring of ECG, heart rate, and respiratory rate.



## Skin Electrode for MEEG Recording

P7001-00296

Skin electrodes for monitoring of MEEG.

25 pcs/package.



## MSpO2 Sensor

P7002-00008

Sensor for registration of pulse and oxygen saturation.

Requires P9001-00501 "MSpO2 Extension Cable" for use.



## MSpO2 Extension Cable

P9001-00501

Extension cable for SpO2 with connector.

Requires P7002-00008 "MSpO2 fingerclip" for use.



## NIBP Cuff Extension Hose

P9001-00485

NIBP extension cable with connector.



## NIBP Cuff Small

P9001-00503

Cuff for measurement of blood pressure (perimeter 20.5-28 cm). Hose length 0.2m.



---

**NIBP Cuff Adult****P9001-00504**

Cuff for measurement of blood pressure (perimeter 27-35 cm). Hose length 0.2m.



---

**NIBP Cuff Large Adult****P9001-00505**

Cuff for measurement of blood pressure (perimeter 34-43 cm). Hose length 0.2m.



---

**NIBP Cuff Adult Thigh****P9001-00506**

Cuff for measurement of blood pressure (perimeter 42-54 cm). Hose length 0.2m.



---

## Consumables

---

**Legplate Belt****CNS000106**

Stretch belt for the legplate.

- 3 cm wide
- Roll of 15 m
- Not made with natural rubber latex



---

**Transducer Belt – Regular****CNS000107**

Stretch belt for the TOCO and FHR transducers.

- 5 cm wide
- Roll of 15 m
- Not made with natural rubber latex



---

**Transducer Belt – Wide****CNS000108**

Stretch belt for the TOCO and FHR transducers.

- 10 cm wide
- Roll of 15 m
- Not made with natural rubber latex



---

**Transducer Belt – Precut****P2224-08001**

Reusable stretch belt for the TOCO and FHR transducers.

- 6 cm wide
- 150 cm long
- Pack of 2 pcs
- Not made with natural rubber latex




---

**Printer Paper (CTG and ST Analysis)****P8105-00063**

Printer paper with CTG+ST analysis grid, with 50-210 bpm FHR scale and 20 bpm per cm scaling.

3 pcs/unit.



---

**Printer Paper (CTG Only)****P8105-00003**

Printer paper with CTG-only grid, with 50-210 bpm FHR scale and 20 bpm per cm scaling.

3 pcs/unit.

